

TV 10

A PIONEER IN THE TREATMENT OF TISSUE LESIONS

TV 10 benefits from 20 years of experience, practice and medical research for an optimal and secure application of vacuumobilization technique.

- TV 10 is the easiest to use and the most complete on technical grounds
- With an extremely simplified conception, TV 10 is reliable and easy to maintain.
- Its unique range of variable depressions -from 0.05 to 1 bar- allows one to treat recent lesions (burn lesions) or even reduce wrinkles or delicate skins while having an anti-fibrosis effect. It does not require the use of tights or any substance which sometimes provoke infections.
- TV 10 has the lowest suction intensity and thus offers the gentlest treatments.
- Research on the techniques of treatment has exhaustively been carried out on a wide range of pathologies with never equalled efficacy and reliability.
- An innovator conception: transparent suction disks allow one to observe the affected area during the treatment process so as to avoid lesions. Their round brims prevents irritation.

The diversity of the diameters has been selected according to the various pathologies to be treated. They include :

- treatment of scars, wrinkles, cellulite
- plastic & reparative surgery
- reflex massages
- sports physiotherapy
- rheumatology

Consequently, they are the most adapted, to each treatment.

- Finally, each user benefits from Eureduc's techniques and experience through an exhaustive training session on how to treat the above pathologies.